



# GeticoFect™ IC Insect Cell Transfection Reagent

## Instruction Manual - Suspension Cells

### Product Overview

GeticoFect™ IC transfection reagent is a proprietary cationic formulation specifically designed for DNA transfection in insect cells. The reagent has been optimized to efficiently transfect various insect cells in conjunction with the baculovirus expression system.

### Composition and Storage

Product Name	Specification	Storage Condition
GeticoFect™ IC Insect Cell Transfection Reagent	0.75 mL / 1.5 mL / 15 mL	Store at 4°C; do not freeze

### Required Materials

#### 1. Provided by This Product

GeticoFect™ IC Insect Cell Transfection Reagent

#### 2. Materials to Be Prepared by the User

Sf9/Sf21 and other insect cells

---

Insect cell culture medium

---

Opti-MEM I reduced-serum medium

---

Bacmid DNA

---

Sterile PETG flat-bottom triangular flask (125 mL)

---



24-well deep well plate

---

27°C CO<sub>2</sub>-free incubator (non-humidified)

---

Orbital shaker platform (for suspension transfection)

---

Optional: Sf9/Sf21 cells

---

#### Other Material List

- Optional: Sf-900 SFM medium
- Optional: 6-well cell culture plate

## Operation Guide

### 1. General Precautions

1. DNA preparation: Purify Bacmid DNA using a dedicated plasmid extraction kit and sterilize it through a 0.22 µm filter before use.
2. Cell resuscitation: Newly thawed cells should be passaged at least twice before use.
3. Reagent mixing: Gently invert the reagent bottle 5–10 times before use to ensure thorough mixing.
4. Complex preparation: Perform in serum-free medium (e.g., Opti-MEM I).
5. Antibiotic contraindication: Do not add antibiotics to the medium during transfection.
6. Cell status: Use healthy cells in the logarithmic growth phase with viability ≥90%.

### 2. Experimental Steps Time Estimation

Step	Time Consumption
Cell preparation	0–60 minutes
Dilution and incubation of transfection reagent	5 minutes
Formation of DNA-transfection reagent complex	5 minutes
Post-transfection culture	72–96 hours

# Transfection Protocols

## 1. Suspension Cell Transfection Protocol (125 mL Shaker Flask)

### 1. Cell Preparation

- Check cell status: Observe cell morphology and density under a microscope to ensure cells are in the logarithmic growth phase with viability  $\geq 90\%$ . Dilute cells to  $(2.5 \times 10^6)$  cells/mL in a total volume of 25 mL using a 125 mL culture flask. Determine cell viability by trypan blue staining: mix cell suspension with equal-volume 0.4% trypan blue, count under a microscope (live cells unstained, dead cells blue).
- Collect cells: Transfer cells (total  $(62.5 \times 10^6)$ ) to a 50 mL centrifuge tube and centrifuge at  $300\times g$  for 5 minutes, placing tubes symmetrically for balance.
- Discard supernatant: Tilt the tube and aspirate supernatant slowly, leaving 0.5–1 mL to protect the cell pellet.
- Resuspend cells: Add 25 mL fresh medium and pipette gently to form a uniform suspension, avoiding bubbles.
- Transfer cells: Move the suspension to a 125 mL ventilated flask, pre-culture on a  $27^\circ\text{C}$  orbital shaker at  $125\pm 5$  rpm (19–25 mm amplitude) or  $95\pm 5$  rpm (50 mm amplitude) for 10–30 minutes.

### 2. Reagent Preparation

- Prepare a sterile 1.5 mL centrifuge tube labeled "GeticoFect™ IC Reagent Dilution".
- Aspirate 30  $\mu\text{L}$  of reagent with a sterile pipette, ensuring correct volume setting.
- Add the reagent to 1 mL of Opti-MEM I medium, invert 5–10 times to mix, and incubate at room temperature for 5 minutes.

### 3. Complex Preparation

- Add 12.5  $\mu\text{g}$  of Bacmid DNA to the diluted reagent, calculate volume based on concentration, invert to mix, and incubate at room temperature for 5 minutes.

### 4. Transfection Operation

- Add the complex dropwise to the flask while gently rotating, then shake the flask to mix cells and complex.
- Culture on the shaker at  $27^\circ\text{C}$  for 72–96 hours, observing cell status (e.g., swelling, suspension as signs of viral infection).



## 2. Adherent Cell Transfection Protocol (6-Well Plate)

### 1. Cell Preparation

- Prepare a cell suspension, count cells to ensure logarithmic phase and viability  $\geq 90\%$ .
- Inoculate  $(1 \times 10^6)$  cells per well in a 6-well plate with 3 mL medium, culture at 27°C for 30–60 minutes for adhesion.
- Check cell adhesion under a microscope before proceeding.

### 2. Reagent Preparation

- Prepare a sterile 1.5 mL tube, aspirate 10  $\mu\text{L}$  of reagent, add to 250  $\mu\text{L}$  Opti-MEM I, invert to mix, and incubate for 5 minutes.

### 3. Complex Preparation

- Add 1  $\mu\text{g}$  of Bacmid DNA to the diluted reagent, mix by inversion, and incubate for 5 minutes.

### 4. Transfection Operation

- Add the complex dropwise to each well, shake the plate gently, and culture at 27°C for 72–96 hours, monitoring cell status.

## Suspension Cell Transfection Scale Adjustment Guide

Parameter	Culture Type			
	Suspension Cells			
Container Type	24-Well Deep Plate	125 mL Flask	250 mL Flask	500 mL Flask
Cell Count Required	$10 \times 10^6$ cells	$62.5 \times 10^6$ cells	$125 \times 10^6$ cells	$250 \times 10^6$ cells
Transfection Culture Volume	4 mL/well	25 mL	50 mL	100 mL
Shaker Speed	250 $\pm$ 5 rpm (19 mm diameter)	125 $\pm$ 5 rpm (19 mm diameter) / 125 $\pm$ 5 rpm (25 mm diameter) / 95 $\pm$ 5 rpm (50 mm diameter)		
Bacmid DNA	1 $\mu\text{g}$	12.5 $\mu\text{g}$	25 $\mu\text{g}$	50 $\mu\text{g}$



Parameter	Culture Type			
GeticoFect™ IC Reagent Volume	5 µL	30 µL	60 µL	120 µL
Opti-MEM Medium Volume	250 µL	1 mL	2 mL	4 mL
Bacmid DNA Volume	2–4 µL	25–50 µL	50–100 µL	100–200 µL

## Notes

- Shaker speed optimization: Determine the optimal speed through pre-experiments as equipment may vary.
- Cytotoxicity: The reagent has low cytotoxicity; medium replacement is unnecessary, but partial replacement can be considered if cell status is poor.
- Storage: Store at 4°C, avoid freezing, and equilibrate to room temperature before use.